



ORTHOPAEDIC
JOINT REPLACEMENT

AESCULAP[®] Centrament[®]

HIP SYSTEM

AESCULAP® Centrament®

PRODUCT CHARACTERISTICS



FEATURES

CONCEPT

There are a wide range of requirements for cemented hip systems today.

With the Centrament stem concept, an implant has been developed that compliments both the improved cementing techniques and the accepted long term results of different cemented prostheses.

DESIGN

The Centrament stem is of a rounded design, this ensures a cohesive cement mantle. The distal guide (centralizer), and the specially designed lateral profile assist with the accurate centering of the Centrament prosthesis stem in the cement bed. Centrament stems are produced from forged CoCr-alloy ISODUR.

SURGERY

The Centrament system is complemented by modular instrumentation, that helps to to achieve a closed cement mantle around the implant.

Combined with the Centrament implants are a range of modular heads and acetabular components, enabling the system to cover a wide range of indications: Total cemented hip replacement, bipolar and hybrid hip arthroplasty.

PRODUCT CHARACTERISTICS



CEMENTED



BIPOLAR



HYBRID



AESCULAP® Centrament®

MODULAR INSTRUMENTATION



Centrament
reamer
cement



Centrament
rasp




Centrament
stem



bone


Lateral wing profile to support stem introduction and load transfer



Rounded surface area reduces cement loading



Distal stem-cross-section provides high rotational stability within the cement mantle



Enhanced cement mantle in the area of load transfer

Smooth prothesis surface promotes uniform cement contact

Conical stem design with proximal increased cement mantle thickness


Narrow cross-section centralizer minimises cement disruption

Centrament = CENTRALIZED IMPLANT POSITION

AESCULAP® Centrament®

SURGICAL PROCEDURE

For preoperative planning, the likely size of the Centrament prosthesis stem and the size of the centralizer should be determined using X-ray templates. In addition to the contour of the prosthesis stem, the templates also contain the outlines of the cement sheath required for anchoring. The outlines correspond to the rasps which must be used to prepare the prosthesis bed.



PREOPERATIVE PLANNING

The standard resection plane is 58° to the shaft axis. As an exception, the resection plane of the Centrament stem 6S is 45° , because it is designed for deformation in dysplasia cases. For intraoperative orientation, implants and instruments are provided with appropriate marks which must conform to the resection plane according to the preoperative plan. The templates contain scales for orientation in the region of the greater trochanter and for planning the resection with orientation relative to the minor trochanter.

Centrament IMPLANTATION



SURGICAL PROCEDURE

Centrament REAMER

The conical reamers are used for preparation of the distal medullary space and are used in ascending order. The reamer 8 is designed for narrow medullary spaces while reamers 10-16 are designed for average ones. The nominal diameter of the largest reamers used corresponds to the distal centralizer to be used.

Centrament CENTRALIZERS

The centralizers consist of PMMA and fit all Centrament stems. The standard size chosen corresponds to the last reamer used. In cases of larger distal medullary space the centralizer can be chosen 2 mm larger.

Centrament RASPS

The rasps are teathed only in their upper part and are used for proximal preparation of the cement-implant-bed. The rasp is centred in the medullary space via the smooth distal part. Trial heads for trial positioning can be mounted on the modular rasp connector.

Centrament STEM SELECTION

The selection of the Centrament stem depends on the last rasp used. The Centrament stems are designed so that, with smaller nominal dimensions than that of the rasp, a cohesive cement mantle is always ensured. The minimum cement thickness on the tip of the prosthesis stem is equivalent to half the difference in nominal diameter from the last rasp used. The cement mantle increases progressively in proximal direction.

AESCULAP® Centrament®

INSTRUMENTS AND IMPLANTS

narrow medullary canal	normal medullary canal					long stem
Ø 8 mm	Ø 10 mm	Ø 12 mm	Ø 14 mm	Ø 16 mm	Ø 14 mm	
 NF491R 8S	 NF492R 10	 NF493R 12	 NF494R 14	 NF495R 16	 NF494R 14	
 NF481R 6S	 NF482R 8	 NF483R 10	 NF484R 12	 NF485R 14	 NF484R 12L	
 NK081K	 NK082K	 NK083K	 NK084K	 NK086K	 NK085K	
 NK088	 NK090	 NK092	 NK094	 NK096	 NK094	
Ø 8 mm	Ø 10 mm	Ø 12 mm	Ø 14 mm	Ø 16 mm	Ø 14 mm	

Centrament stems



ISODUR_F

6S	NK081K	140 mm
8	NK082K	150 mm
10	NK083K	155 mm
12	NK084K	160 mm
12L	NK085K	220 mm
14	NK086K	165 mm

Centralizers



PMMA

8 mm	NK088
10 mm	NK090
12 mm	NK092
14 mm	NK094
16 mm	NK096

Modular heads



12/14
ISODUR_F

Ø	22.2 mm	28 mm	32 mm
S	–	NK429K	NK529K
M	NK330K	NK430K	NK530K
L	NK331K	NK431K	NK531K
XL	–	NK432K	NK532K
XXL	–	NK433K	NK533K

Implant materials:

- ISODUR_F
Cobalt-chromium forged alloy
(CoCrMo/ISO 5832-12)
- BioloX delta
Aluminiumoxyd–Matrix–Ceramic
(Al₂O₃/ZrO₂/ISO 6474-2)
- PMMA
Polymethylmethacrylate



12/14
BioloX delta

Ø	22.2 mm	28 mm	32 mm
S	–	NK460D	NK560D
M	–	NK461D	NK561D
L	–	NK462D	NK562D
XL	–	–	NK563D

BioloX delta

AESCULAP® Centrament®

INSTRUMENT SETS



NF 500 Centrament rasps consisting of:

1	Centrament rasp 8S	NF481R
1	Centrament rasp 10	NF482R
1	Centrament rasp 12	NF483R
1	Centrament rasp 14	NF484R
1	Centrament rasp 16	NF485R
2	Modular rasp handle	NG115R
1	Tray for form rasps	NF499R
1	Wrapping cloth	JF511
1	Identification plate	JG645B

Please order separately:

1	Trial prosthesis head 12/14 22.2 mm M	NF327
1	Trial prosthesis head 12/14 22.2 mm L	NF328
1	Trial prosthesis head 12/14 28 mm S	NF336
1	Trial prosthesis head 12/14 28 mm M	NF337
1	Trial prosthesis head 12/14 28 mm L	NF338
1	Trial prosthesis head 12/14 28 mm XL	NF339
1	Trial prosthesis head 12/14 28 mm XXL	NF343
1	Trial prosthesis head 12/14 32 mm S	NF346
1	Trial prosthesis head 12/14 32 mm M	NF347
1	Trial prosthesis head 12/14 32 mm L	NF348
1	Trial prosthesis head 12/14 32 mm XL	NF349
1	Trial prosthesis head 12/14 32 mm XXL	NF353

Recommended container for NF500 and NF502 Aesculap basic container 592 x 285 x 153 mm

X-ray-templates (please order separately):

Centrament size 8+10	NF081
Centrament size 12+14	NF082
Centrament 6S	NF083
Centrament 12L	NF085



NF 502 Centrament instruments consisting of:

1	Centrament form reamer size 8S Harris	NF491R
1	Centrament form reamer size 10 Harris	NF492R
1	Centrament form reamer size 12 Harris	NF493R
1	Centrament form reamer size 14 Harris	NF494R
1	Centrament form reamer size 16 Harris	NF495R
1	T-handle L 125 mm w/Harris chuck	ND144R
1	Preparation reamer D 8.0 mm w/T-handle	ND359R
1	Stem impactor	ND830R
1	Centrament insertion handle	ND824R
1	Centrament tray instruments/ form reamers	NF501R
1	Wrapping cloth	JF511
1	Identification plate	JG645B

Please order separately:

1	Insertion instrument for medullary bone plugs	NG702R
1	Intermedullary bone plug trephine diam. 8-10 mm Harris	ND185R
1	Intermedullary bone plug trephine diam. 10-12.5 mm Harris	ND186R
1	Intermedullary bone plug trephine diam. 12.5-15 mm Harris	ND187R
1	Intermedullary bone plug trephine diam. 15-18 mm Harris	ND189R

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" and the product trademarks "Centrament" and "ISODUR" are registered trademarks of Aesculap AG. "BioloX" is a registered trademark of CeramTec GmbH, Plochingen.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.