


# Aesculap<sup>®</sup> Plasmacup<sup>®</sup>

Cementless Acetabular Cup System




Aesculap Orthopaedics



Plasmacup® SC



Plasmacup® delta



# Content

■ Plasmacup® System	5
■ Plasmapore® Surface	6
■ Plasmacup® Design	7
■ Plasmacup® Articulation	8
■ BioloX® Ceramic-on-Ceramic	9
■ BioloX® delta Ceramics	10
■ Plasmacup® delta Concept	12
■ Plasmacup® Surgical Technique	14
■ OrthoPilot® THA Navigation	19
■ Plasmacup® Ordering Information	20
Instruments	22
Implants	26



# Plasmacup® System

The Plasmacup® system, which has proven itself successful in clinical practice since 1992, based on three essential elements:




Plasmapore®

- Microporous Plasmapore® titanium coating for excellent primary and secondary stability




Plasmacup®

- Successful cup system with ceramic-ceramic articulation



OrthoPilot®

- One of the world's leading navigation technology in hip arthroplasty



BioloX® delta

- Low wear and high demand treatment in total hip arthroplasty


# Plasmapore® Surface



Plasmapore®



Primary stability



Plasmapore® μ-CaP

Plasmapore® coated implants have been used by Aesculap since 1986. In a vacuum coating process, pure titanium powder is applied to the surface of cementless implants, to form a 0.35 mm thick layer with up to 40 % microporosity.

The pore size of the Plasmapore® coating ranges between 50 and 200 μm to allow direct bone apposition.


The rough surface of the Plasmapore® structure supports a more stable primary fixation of the implants than other implant surface types. Stability measurements demonstrate the superior primary stability achieved with Plasmapore® compared with other implant surfaces.

The newly developed Plasmapore® μ-CaP coating uses calcium phosphate as a coating material.

A 20 μm layer of high-purity dicalcium phosphate dihydrate (DCPD) is electrochemically applied to the Plasmapore® coating. The thin μ-CaP coating accelerates the formation of bone material at the implant surface and dissolves, within 8 to 12 weeks, without the involvement of macrophages.

The Plasmacup® SC implants are available with Plasmapore® or Plasmapore® μ-CaP surfaces.

Further detailed information on Plasmapore® μ-CaP can be found in the Aesculap brochure 051002.



Primary stability of various press-fit cups in biomechanical experiments

Pitto RP, Bohner J, Hofmeister V.  
Factors affecting the primary stability of acetabular components.  
An in vitro study *Biomed Tech (Berl)*. 1997 Dec;42(12):363-8.


# Plasmacup® Design



Design



Stability



Inner surface


The Plasmacup® is characterized by good press-fit stability and safe attachment of modular polyethylene or ceramic liners.

The external shape of Plasmacup® is hemispherical with a slightly flattened apex. With standard liners the center of rotation is located at the exact center of the sphere. Plasmacup® is suitable for either polyethylene (PE) or ceramic cup liners.

The Plasmacup® liners are attached by a press-fit cone with a large surface area and, in case of the PE liners, throughfull contact with the base of the cup. In this way, both polyethylene and ceramic liners are well fixed. The drill holes are located in the cranial region of the cup, outside the conical attachment surface. The rough titanium inner surface reduces relative movements to only a few microns, which prevents the formation of abrasion particles on the back side of the liner.

The conical fixation surface of the Plasmacup® polyethylene liners also forms a seal against the migration of polyethylene particles from the articulating joint, and thus reduces the risk of osteolysis adjacent to the screw holes.

The polyethylene liners are strongest when the load is directed cranially. In the primary load area, Plasmacup® polyethylene implants are at least 6 mm thick. The fixation is highly stable against tilting and rotation forces in vivo.



Overview of the most important data regarding fixation, for Plasmacup® PE and BioloX® liners

Blömer W.  
Design aspects of modular inlay fixation  
Hip International 1997, Vol. 7, No. 3:110-20.

# Plasmacup® Articulation



System SC



Polyethylene



BioloX®


Plasmacup® implants can be implanted with polyethylene or BioloX® ceramic liners. Implants with ceramic liners are marked with the letters SC (System Ceramic).

The polyethylene used by Aesculap conforms to established standards and long-term clinical experience. The implants are manufactured from high density PE plates, using CNC technology. The material is sterilized by radiation in a nitrogen atmosphere, a process that has been established in Europe since the mid-eighties.

Modern packaging materials protect the polyethylene implants against oxygen during and after sterilization. In-vivo wear of the Plasmacup® polyethylene liners with a ceramic 28 mm head is 0.1 mm per year, which is below the threshold that would cause osteolysis in terms of the number of the PE particles. Higher wear can occur with metal heads, by third-body wear, through incorrect cup positioning or as a result of implant loosening.

Implantation of ceramic Plasmacup® liners made with BioloX® reduces wear in the joint to a few µm per year. This implant, correctly positioned and with stable fixation, is widely used in the treatment of young patients.

Plasmacup® is one of the leading implant systems with BioloX® ceramic liners.




ISO 14242 hip simulator wear measurements and data referring to other studies

Graphic/Source: Dr. Ing. Christian Kaddick, Endolab Mechanical Engineering GmbH, Thansau/Rosenheim.



# Bilox<sup>®</sup> Ceramic-on-Ceramic



Plasmacup<sup>®</sup> SC was specially designed for the use of ceramic Bilox<sup>®</sup> liners. The conical inner shape of Plasmacup<sup>®</sup> is also used for attaching polyethylene liners, which means that the surgeon is always free to choose the bearing material that is best for the patient.


28 mm Bilox<sup>®</sup> liners are available for Plasmacup<sup>®</sup> implants from size 44 mm, 32 mm Bilox<sup>®</sup> liners for implants from size 48 mm and 36 mm Bilox<sup>®</sup> liners for Plasmacup<sup>®</sup> SC implants from size 56 mm.

Range of motion and dislocation stability of the hip implant depend on the head diameter and the trunnion size of the prosthesis. For the ceramic-ceramic bearing surfaces, Aesculap recommends using 32 mm and 36 mm heads. Additionally, stems with 8/10 trunnion are available, which enhances the implant range of motion even more, up to 150 degrees with a 36 mm head.

Since the design of the ceramic-ceramic THA articulation does not provide anti-dislocation elements, any tendency towards subluxation or dislocation of the joint

constitutes a contraindication for this articulating surface.

The implantation of ceramic liners is also contraindicated in cases of a socket position of more than 55° inclination, retroversion or excessive anteversion of more than 25°. Such implant situations can lead to excessive load on the articular surface of the prosthesis head and on the rim of the cup liner. This load cannot be corrected, even through implanting a larger head.



*Hip range of motion with different head diameters and prosthesis cone sizes*




Bilox<sup>®</sup> implants


# Bilox® delta Ceramics



BioloX® delta





BioloX® forte



BioloX® delta

BioloX® delta, the latest generation of ceramics in THA, is a high strength aluminium oxide matrix ceramic. Finest  $ZrO_2$  particles represent the matrix material. This leads to an increased material strength in direct comparison with BioloX® forte. BioloX® delta prosthesis heads and inserts can articulate with each other and without any limitations also with BioloX® forte implant components. BioloX® delta offers a higher implant strength while keeping the excellent ceramic wear characteristics. Therefore cup inserts for larger head diameters and prosthesis heads with XL neck length can be provided with BioloX® delta ceramics. The BioloX® delta ceramic articulation sets a new benchmark for high demand hip arthroplasty.



## BioloX® delta 36 mm Ceramic-on-Ceramic THA



BioloX® delta is an innovative material development for a low wear high demand hip joint replacement.

BioloX® delta reduces the risk of any articulation failure. High strength, larger head diameters and new developed acetabular components will achieve and contribute to a new generation of ceramic implant components.

The new Plasmacup® delta implants extend the ratio of surgeries which are preoperatively planned and indicated for a 36 mm ceramic on ceramic THA for cup sizes 52 and 54.

Plasmacup® delta and 36 mm ceramic BioloX® delta insert are delivered in one package. For cup sizes 52 and 54 as a modular system and for smaller cup sizes 48 and 50 preassembled.

Plasmacup® delta features a thinner shell thickness. The Plasmapore® coating covers the complete cup surface as there are no additional holes for screw fixation. Plasmacup® delta is implanted with the instruments of the Plasmacup® system using the same surgical technique. The preassembled implants are inserted with special cup attachments. The implantation can be naturally navigated with OrthoPilot®.



# Plasmacup® delta Concept


The future of ceramic THA articulation

- Large 36 mm head diameter
- Biolox® delta high strength Al<sub>2</sub>O<sub>3</sub>-matrix-composite-ceramic
- Plasmacup® delta preassembled for smaller cup sizes
- Plasmacup® delta modular for cup sizes 52 and 54
- Extra-long neck for ceramic heads






Plasmacup® delta 44 mm



Plasmacup® delta 48 mm



Plasmacup® delta 52 mm

Plasmacup® delta enables a large ceramic articulation for small cup diameters and extends therefore the indication range of Plasmacup® SC.

Plasmacup® delta is suitable for larger head diameters for improved joint stability and increased range of motion.

Plasmacup® delta has a thinner metal shell compared to Plasmacup® SC with a closed Plasmapore® surface without screw holes.

For cup sizes 52 and 54 a modular Plasmacup® delta is implanted with the standard Plasmacup® system instruments using the same surgical technique.

The ceramic liner is packaged together with the Plasmacup® delta shell and is inserted after shell implantation and trial reduction.

The small, semi-modular implant solutions for cup sizes 44 – 50 mm are inserted with special cup attachments. The implants are delivered preassembled. In revision cases the ceramic inserts can be removed.

The implantation procedure can be navigated with OrthoPilot®.



# Plasmacup® Surgical Technique



Acetabular reamers



Trial cups



Plasmacup®

The Plasmacup® implantation instruments have undergone continuous development during more than 10 years of clinical application. The optional use of the OrthoPilot® hip navigation system sets the trend for a safe and reproducible operating technique (see p. 19).

In order to have a good pressfit fixation of the Plasmacup® implant, there must be a good bony structure and proper surgical technique. Acetabular exposure removal of the articular cartilage and osteophytes are required for the proper preparation of the acetabulum.

This is done using spherical reamers, which are driven by a low-speed motor handpiece. During the reaming procedure, all cartilaginous material must be ablated down to the subchondral bone until bleeding occurs.

For non-dysplastic cases, care must be taken that the center of rotation of the joint is not medialized unnecessarily. The socket edges should be prepared for a sufficiently large bony fixation surface.


In cases of dysplastic changes, a cup position in the region of the primary socket is recommended, as far as a shortening of

the leg can be compensated. The caudal edge of the socket should be at the level of the tear drop figure. A cranial bone graft is performed, if necessary, before the socket base is deepened to provide sufficient cranial roofing.

The size of the Plasmacup® implant corresponds to the size of the last acetabular reamer used and includes the proper pressfit conditions.

The final selection of the implant is only determined after a trial cup has been seated firmly. A stable fit of this trial cup


is achieved when the pelvis of the patient can be moved by gently moving the cup impactor by about 10 degrees. The trial implant can be easily levered out from the in-vivo trial position by moving beyond this angle.



Reaming of the acetabulum


Plasmacup® trial cup

## Plasmacup® Cup Position




*Aiming device supine position*

Aiming devices are available for Plasmacup® to measure inclination and anteversion for both standard and navigated surgeries. These devices, which have been designed for supine or lateral decubitus position, can be mounted on the cup impactor shaft.



*Aiming device lateral decubitus position*


The accurate and stable assembly of the Plasmacup® implant on the impactor shaft must be checked by a surgical assistant and the surgeon prior to implantation. A slotted hammer on the impactor shaft is suitable for shifting and correcting the position of the Plasmacup® implant.



*Plasmacup® with trial liner*

After completing the Plasmacup® surgery steps of acetabular exposure, reaming, assessment of the cup bed with the trial implant and implantation of the cup (Plasmacup® SC, NSC or MSC), the trial liner is inserted. The final selection of the modular liner (PE or ceramics) is determined only after the stem is implanted and a final trial reduction has been performed.

# Plasmacup® Implantation



*Plasmacup® fixation screws*

*The pivot angle of the 6.5 mm Plasmacup® screws is 20 degrees. Before implanting the modular liner the surgeon has to make sure that none of the screw heads protrude into the liner anchoring zone.*

Generally, in good bone Plasmacup® can be implanted without additional screws. As a stability check the cup impactor is moved through  $\pm 20$  degrees until the patient's pelvis moves. Under these conditions, Plasmacup® SC can also be rotated 180° prior to implantation, placing the screw holes in the non load bearing caudal region since they are not needed cranially. If there is any doubt concerning the intraoperative primary stability, fixation screws must be used or the implant must be replaced.

For cases where additional stability with fixation screws is necessary, Plasmacup® SC features three holes in the cranial region. To protect the medial blood vessels, the middle and lateral screw positions can be used and the medial hole is usually left open. Plasmacup® MSC offers further additional screw holes in the cranial and caudal region. The cranial holes are located further laterally or medially with this implant.

Prior to inserting the self-tapping 6.5 mm screws, the drill holes are prepared with a flexible 3.2 mm drill (or with a 4 mm drill for severely sclerotic bone). The required screw length is measured and the screws are implanted using a screw holding forceps and a cardan-jointed screwdriver.



*Instruments for screw implantation*



# Plasmacup® Liners



PE standard liner




PE posterior wall liner



PE asymmetrical 10 deg. liner



BioloX® delta ceramic liner



*Fingertip check of ceramic liner position*

Posterior wall (hooded) Plasmacup® PE liners increase luxation stability e.g. towards posterior for implantations using the posterior surgical approach. The asymmetrical liners correct the cup position by 10 degrees.

When using ceramic liners the final check for seating is assessed with a fingertip check. After seating, the liner is fixed using an impactor with a plastic head.



The ceramic Plasmacup® liners can be removed with a punch. When doing this it is important to place the instruments correctly in the pits on the edge of the implant, and separate the liner from the cup with several sharp blows or impulses\*.

*\* Also see the instructions for use enclosed with every Plasmacup® implant.*

*Instruments for removing ceramic liners*

# Plasmacup® delta Surgical Technique



*Prepare insertion instrument*

For the implantation of the preassembled Plasmacup® delta cup implants sizes 44 – 50 a special insertion attachment has to be used. For the sizes 44 and 46 mm a device with 32 mm diameter has been created (ND266R), for the sizes 48 and 50 mm an attachment with 36 mm diameter (ND270R). The additional insertion device is mounted on the standard Plasmacup® impactor. The rim of the attachment pushes during the impaction process on the titanium shell of the Plasmacup® delta and prevents at the same time the loosening of the premounted ceramic insert.



*Fill in water*

The connection between semi-modular (preassembled) Plasmacup® delta and the impactor is realized by vacuum. Semi-modular indicates the possibility to remove the ceramic insert in any case of revision. To create the vacuum the ceramic insert has to be filled completely with water. The impactor with the attachment is put on the implant. With pushing in and twisting at the same time the plastic part of the attachment device sucks on the wet ceramic insert.



*Push in and twist*

During the impacting process the implant can only be guided in implanting direction. With light pressure on the instrument, and thus on the rim of the Plasmacup® delta, the implanting direction can be influenced slightly. An extraction is not possible because the suction power is not strong enough. Therefore it is mandatory to check the prepared acetabular implant bed with the equivalent Plasmacup® trial cup before implantation.

*Note: The insertion attachment can be used with OrthoPilot® navigation*

# OrthoPilot® THA Navigation



OrthoPilot® transmitter referencing

All Plasmacup® components can be used with OrthoPilot® navigation technology. OrthoPilot® cup navigation works without CT or fluoroscopy, following the principles of kinematic referencing.

In navigated Plasmacup® surgeries the system measures the inclination and anteversion angles relative to the anterior pelvic plane. During the acetabular reaming stage, the joint center, the reaming depth and the orientation of the reamer are measured and displayed.



OrthoPilot® navigation of acetabular reaming

Plasmacup® navigation with OrthoPilot® is suitable for different patient positions and surgical approaches. It also supports less invasive surgical procedures and surgeries on dysplastic cases. Navigation of the cup is an integral part of OrthoPilot® THA navigation.




OrthoPilot® hip navigation of leg length and offset

The new OrthoPilot® THA Plus navigation combines the data on the position of Plasmacup® with the position of the stem. The surgeon obtains information regarding leg length and offset of the hip joint. Especially in less invasive procedures OrthoPilot® THA Plus supports the surgeon during surgery. Further information on OrthoPilot® hip navigation can be found in the hip navigation brochure O21902.



# Plasmacup® Acetabulum Reamers



## Basket tray NF932R

Aesculap basket tray 485 x 253 x 76 mm with supports for:

- 13 reamers (e.g. 44 to 68 mm)
- 2 straight reamer shanks (e.g. FS960R)
- 1 curved reamer shank (e.g. NF936R)
- Standard sleeve FS974 and navigation sleeve FS939

NF933R Aesculap basket tray 485 x 253 x 76 mm (w/o fig.) supports for:

- 24 reamers
- 2 straight reamer shanks (e.g. FS960R)



Please order separately					
∅		∅		∅	
40 mm	NF940R	50 mm	NF950R	60 mm	NF960R
42 mm	NF942R	52 mm	NF952R	62 mm	NF962R
44 mm	NF944R	54 mm	NF954R	64 mm	NF964R
46 mm	NF946R	56 mm	NF956R	66 mm	NF966R
48 mm	NF948R	58 mm	NF958R	68 mm	NF968R



Please order separately	
OrthoPilot® reamer shank ZIMMER	FS959R
OrthoPilot® reamer shank Harris	FS960R
OrthoPilot® reamer shank AO	FS961R
OrthoPilot® sleeve for FS959R to FS961R	FS939
Standard sleeve for FS959R to FS961R	FS974

Please order separately	
Curved reamer shank ZIMMER	NF935R
Curved reamer shank Harris	NF936R
Curved reamer shank AO	NF937R

# Aesculap Classic Reamers




Acetabulum reamers			
40 mm	NG540R	56 mm	NG556R
42 mm	NG542R	58 mm	NG558R
44 mm	NG544R	60 mm	NG560R
46 mm	NG546R	62 mm	NG562R
48 mm	NG548R	64 mm	NG564R
50 mm	NG550R	66 mm	NG566R
52 mm	NG552R	68 mm	NG568R
54 mm	NG554R		

Reamer shanks	∅ 40 - 48 mm	∅ 50 - 68 mm
Harris	NG621R	NG631R
AO	NG623R	NG633R
triangular	NG627R	NG637R
Hudson	NG629R	NG639R

Note: The Aesculap classic reamers are only delivered as replacements

# Plasmacup® Instruments



## Plasmacup® instrument set 1 NF240

Comprising		Please order separately	
Insertion instrument straight	FS944R	1/1 size wide perforated basket lid	JH217R
Tray for NF240 (48 x 253 x 76 mm)	NF241R	Aiming device for supine position	NF277R
Grafic template for NF241R (NF240)	TE912	Aiming device for posterior approach	NF292R
Cloth for lining deep containers	JF511	Insertion instrument curved	FS947R
		T-Handle for insertion instrument	FS948R
		Screw driver for FS947R	NF371R
		Cup pressing head ø 22.2 mm	ND178
		Cup pressing head ø 26 mm	ND179
		Cup pressing head ø 28 mm	ND174
		Cup pressing head ø 32 mm	ND172
		Cup pressing head ø 36 mm	ND166




#### Plasmacup® instrument set 2 NF242

Comprising	
Cup inserting and pressing instrument	ND170R
Impaction/extraction instrument	ND401R
Slotted hammer	NF275R
Removal forceps for PE-inserts	NG430R
Articulated screw driver SW 3.5	NF285R
Screw holding forceps	NF287R
Screw gauge	NF269R
Drill guide for screw $\varnothing$ 3.2 mm	NF278R
Drill guide for screw $\varnothing$ 4.0 mm	NF279R
Flexible drill $\varnothing$ 3.2/32 mm	NF280R
Flexible drill $\varnothing$ 3.2/44 mm	NF281R
Flexible drill $\varnothing$ 4.0/32 mm	NF282R
Tray for NF242 (48 x 253 x 76 mm)	NF243R
Grafic template for NF243R (NF242)	TE913
Cloth for lining deep containers	JF511

Please order separately	
1/1 size wide perforated basket lid	JH217R

# Plasmacup® Instruments



## Plasmacup® trial cups and trial inserts NG036

### Comprising

#### Plasmacup® SC/MSC trial cup

ø		ø	
44 mm	NG944R	58 mm	NG958R
46 mm	NG946R	60 mm	NG960R
48 mm	NG948R	62 mm	NG962R
50 mm	NG950R	64 mm	NG964R
52 mm	NG952R	66 mm	NG966R
54 mm	NG954R	68 mm	NG968R
56 mm	NG956R		

Removal forceps for trial inserts NG437R

Tray for NG036 (489 x 253 x 48 mm) NG037R

Cloth for lining deep containers JF511

### Please order separately

1/1 size wide perforated baset lid JH217R

Silicone basket liner fitting JF159R JF946

Plasmacup® SC/MSC trial cup size 40 mm NG940R

Plasmacup® SC/MSC trial cup size 42 mm NG942R

Recommended container for NG036 and the acetabulum reamer set Aesculap basic container 592 x 274 x 187 mm (e.g. JK444)



Please order separately

Trial liners	symmetrical				posterior wall			asymmetrical	
	ø 22.2 mm	ø 28 mm	ø 32 mm	ø 36 mm	ø 22.2 mm	ø 28 mm	ø 32 mm	ø 28 mm	ø 32 mm
ø 40/42 mm	NG370	–	–	–	NG600	–	–	–	–
ø 44/46 mm	NG371	NG391	–	–	NG601	NG641	–	NG491	–
ø 48/50 mm	NG372	NG392	NG502	–	NG602	NG642	–	NG492	–
ø 52/54 mm	NG373	NG393	NG503	–	NG603	NG643	NG513	NG493	NG573
ø 56/58 mm	NG374	NG394	NG504	NG509	NG604	NG644	NG514	NG494	NG574
ø 60/62 mm	NG375	NG395	NG505	NG510	NG605	NG645	NG515	NG495	NG575
ø 64 - 68 mm	NG376	NG396	NG506	NG511	–	NG646	NG516	NG496	NG576

# Plasmacup® SC Implants



## Plasmacup® SC

40 mm	NH040T
42 mm	NH042T
44 mm	NH044T
46 mm	NH046T
48 mm	NH048T
50 mm	NH050T
52 mm	NH052T
54 mm	NH054T
56 mm	NH056T
58 mm	NH058T
60 mm	NH060T
62 mm	NH062T
64 mm	NH064T
66 mm	NH066T
68 mm	NH068T

ISOTAN®<sub>F</sub>

## Plasmacup® SC μ-CaP


40 mm	—
42 mm	—
44 mm	NC444T
46 mm	NC446T
48 mm	NC448T
50 mm	NC450T
52 mm	NC452T
54 mm	NC454T
56 mm	NC456T
58 mm	NC458T
60 mm	NC460T
62 mm	NC462T
64 mm	NC464T
66 mm	NC466T
68 mm	NC468T

ISOTAN®<sub>F</sub>

## Plasmacup® NSC

40 mm	NH340T
42 mm	NH342T
44 mm	NH344T
46 mm	NH346T
48 mm	NH348T
50 mm	NH350T
52 mm	NH352T
54 mm	NH354T
56 mm	NH356T
58 mm	NH358T
60 mm	NH360T
62 mm	NH362T
64 mm	NH364T
66 mm	NH366T
68 mm	NH368T

ISOTAN®<sub>F</sub>



### Plasmacup® MSC

40 mm	NH140T
42 mm	NH142T
44 mm	NH144T
46 mm	NH146T
48 mm	NH148T
50 mm	NH150T
52 mm	NH152T
54 mm	NH154T
56 mm	NH156T
58 mm	NH158T
60 mm	NH160T
62 mm	NH162T
64 mm	NH164T
66 mm	NH166T
68 mm	NH168T

ISOTAN®<sub>F</sub>

### Plasmacup® MSC μ-CaP

40 mm	NC540T
42 mm	NC542T
44 mm	NC544T
46 mm	NC546T
48 mm	NC548T
50 mm	NC550T
52 mm	NC552T
54 mm	NC554T
56 mm	NC556T
58 mm	NC558T
60 mm	NC560T
62 mm	NC562T
64 mm	—
66 mm	—
68 mm	—

ISOTAN®<sub>F</sub>

### Plasmacup® screws 6.5 mm

16 mm	NA766T
20 mm	NA770T
24 mm	NA774T
28 mm	NA778T
32 mm	NA782T
36 mm	NA786T
40 mm	NA790T
44 mm	NA794T

ISOTAN®<sub>F</sub>

# Plasmacup® SC Liners



Plasmacup® SC cup liners  
Polyethylene

	symmetrical			posterior wall			asymmetrical	
	ø 22.2 mm	ø 28 mm	ø 32 mm	ø 22.2 mm	ø 28 mm	ø 32 mm	ø 28 mm	ø 32 mm
40 mm 42 mm	NH170	–	–	NH300	–	–	–	–
44 mm 46 mm	NH171	NH191	–	NH301	NH401	–	NH471	–
48 mm 50 mm	NH172	NH192	NH202	NH302	NH402	–	NH472	–
52 mm 54 mm	NH173	NH193	NH203	NH303	NH403	NH413	NH473	NH323
56 mm 58 mm	NH174	NH194	NH204	NH304	NH404	NH414	NH474	NH324
60 mm 62 mm	NH175	NH195	NH205	NH305	NH405	NH415	NH475	NH325
64 mm 66 mm 68 mm	NH176	NH196	NH206	NH306	NH406	NH416	NH476	NH326

UHMWPE

## Plasmacup® X-ray templates

	symmetrical	posterior wall	asymmetrical
44 - 52 mm	NG400	NG418	NG403
54 - 62 mm	NG401	NG419	NG404
64 - 68 mm	NG402	NG420	–



### Ceramic

	ø 28 mm	symmetrical ø 32 mm	ø 36 mm
40 mm 42 mm	–	–	–
44 mm 46 mm	NH091D	–	–
48 mm 50 mm	NH092D	NH102D	–
52 mm 54 mm	NH093D	NH103D	–
56 mm 58 mm	–	NH104D	NH109D
60 mm 62 mm	–	NH105D	NH110D
64 mm 66 mm 68 mm	–	NH106D	NH111D

BioloX® delta

# Plasmacup® delta



## Plasmacup® delta

		liner
44 mm	NH644D	ø 32 mm semi-modular
46 mm	NH646D	ø 32 mm semi-modular
48 mm	NH648D	ø 36 mm semi-modular
50 mm	NH650D	ø 36 mm semi-modular
52 mm	NH652D	ø 36 mm modular
54 mm	NH654D	ø 36 mm modular

Plasmacup® delta implants complete the Plasmacup® SC program with 36 mm ceramic liners. These implants can not be combined with Plasmacup® SC components and are only supplied together with the BioloX® delta cup component. Special liners with shoulder are available for revision operations.

BioloX® delta

## Plasmacup® delta instruments

Please order separately	
44/46 attachment ø 32 mm	ND266R
48/50 attachment ø 36 mm	ND270R
Replacement o-ring for ND266R	TA012710
Replacement o-ring for ND270R	TA012640



**Plasmacup® delta liners for revisions**  
Polyethylene

	posterior wall	
	ø 28 mm	ø 32 mm
44 mm	NH407	–
46 mm		
48 mm	–	NH417
50 mm – 54 mm	–	NH418

UHMWPE



**Ceramic**

	symmetrical	
	ø 32 mm	ø 36 mm
44 mm	NH632D	–
46 mm		
48 mm	–	NH636D
50 mm – 54 mm	–	NH637D

BioloX® delta

**Implant Materials:**

ISOTAN® F Titanium forged alloy (Ti6Al4V / ISO 5832-3)

Plasmapore® Pure titanium (Ti / ISO 5832-2)

Plasmapore® µ-CaP Pure titanium surface with 20 µm coating dicalcium phosphate dihydrate (CaHPO<sub>4</sub> x 2H<sub>2</sub>O)

BioloX® delta Aluminium oxide matrix composite ceramic (Al<sub>2</sub>O<sub>3</sub> / ZrO<sub>2</sub> / ISO 6474-2)

UHMWPE Ultra high molecular weight polyethylene (ISO 5834-2)

