


Excia CEMENTLESS		Size	CCD	Stem Length (A)	Stem Length (F)	Offset (B)	Neck Length [C]	Neck Height (D)	Proximal (G)	Distal (Size) (H)	Height (J)
Standard	8	135°		112	145	38	33	30	28	8	18
	9	135°		116	150	39	34	31	29	9	19
	10	135°		120	154	40	35	32	30	10	20
	11	135°		124	159	41	36	32	31	11	21
	12	135°		127	163	43	38	33	32	12	22
	13	135°		131	168	44	39	34	34	13	23
	14	135°		135	172	45	40	35	35	14	24
	15	135°		138	177	46	41	36	36	15	25
	16	135°		142	181	47	42	36	37	16	26
	17	135°		146	186	49	43	37	38	17	27
18	135°		150	190	50	44	38	39	18	28	
Lateralisiert	8	128°		112	145	44	38	30	28	8	18
	9	128°		116	150	45	39	31	29	9	19
	10	128°		120	154	46	40	32	30	10	20
	11	128°		124	159	47	41	32	31	11	21
	12	128°		127	163	49	42	33	32	12	22
	13	128°		131	168	50	43	34	34	13	23
	14	128°		135	172	51	45	35	35	14	24
	15	128°		138	177	52	46	36	36	15	25
	16	128°		142	181	53	47	36	37	16	26
	17	128°		146	186	55	48	37	38	17	27
18	128°		150	190	56	49	38	39	18	28	
Excia CEMENTED		Size	CCD	Stem Length (A)	Stem Length (F)	Offset (B)	Neck Length [C]	Neck Height (D)	Proximal (G)	Distal (Size) (H)	Height (J)
Standard	8	135°		-	-	38	33	30	-	-	-
	9	135°		107	140	39	34	31	23	7	9
	10	135°		110	144	40	35	32	24	8	10
	11	135°		114	148	41	36	32	25	9	11
	12	135°		118	153	43	38	33	26	10	12
	13	135°		122	158	44	39	34	27	11	13
	14	135°		125	162	45	40	35	28	12	14
	15	135°		129	167	46	41	36	30	13	15
	16	135°		133	171	47	42	36	31	14	16
	17	135°		136	175	49	43	37	32	15	17
18	135°		140	180	50	44	38	33	16	18	
Lateralisiert	8	128°		-	-	44	38	30	-	-	-
	9	128°		107	140	45	39	31	23	7	9
	10	128°		110	144	46	40	32	24	8	10
	11	128°		114	148	47	41	32	25	9	11
	12	128°		118	153	49	42	33	26	10	12
	13	128°		122	158	50	43	34	27	11	13
	14	128°		125	162	51	45	35	28	12	14
	15	128°		129	167	52	46	36	30	13	15
	16	128°		133	171	53	47	36	31	14	16
	17	128°		136	175	55	48	37	32	15	17
18	128°		140	180	56	49	38	33	16	18	



Date of issue	08 10 2020	IMPORTANT REMARKS
Place of issue	Tuttlingen, Germany	Values are reflecting the dimensions of the mentioned stem above Values refer to Aesculap drawings and belonging stem characteristics
Dpt.	Dominik Kuehne Director Global Marketing Hip Athroplasty	CHANGES OF LISTED DESIGN PARAMETERS ARE STRICTLY PROHIBITED

GENEHMIGT
 Von Admin, 14.05, 18/9/2009

GENEHMIGT
 Von Admin, 14.05, 18/9/2009